

## **Small but strong IR-15 Wrist Transmitter**

IR-1S Wristband transmitter is an IR remote control with only one button intended to assist persons who are unable to use a standard remote control due to motor, tactile or visual impairments.

The device can help the user to control electronic equipment at home, such as lighting and door openers. The device is intended for indoor and outdoor use. The bracelet transmitter can be used with the included bracelet or necklace.

IR-1S Wristband transmitter is a 1-channel IR transmitter that can be programmed with 7 different channels.

The product is marketed as a technical aid class I for people with disabilities and complies with Regulation (EU) 2017/745 - MDR.

Art nr 400500



## **Technical information:**

## **Accessories:**

Physical goals: (L x B x H): 46 x 44 x 14 mm Battery pack, art nr 401930

Weight incl batteries: 25 gram

